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**INTRODUCTION**

This document contains the precautions, and instructions on how to use and treat (cleaning, disinfecting and sterilizing) all ancillary devices manufactured by DEDIENNE Santé. These instructions have been validated and allow effective treatment of ancillary devices, so they can be reused.

They have been drafted based on the following standards:

- NF EN ISO 16061: 2015 "Instrumentation for use in association with non-active surgical implants - General requirements"
- NF EN ISO 17664: 2017 "Processing of health care products - Information to be provided by the medical device manufacturer for the processing of medical devices»

The instruments covered in this document are reusable instruments for use during implantation of orthopaedic implants manufactured by DEDIENNE Santé.

**SYMBOLS AND ABBREVIATIONS**

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	"NB, see Instructions for Use" and "Consult the supporting documentation"
	Batch number, ensuring product traceability
	Catalogue reference
	CE mark followed by ID number of the notified body for class IIa ancillary devices
	CE mark for class I ancillary devices
	Device manufacturer, with name and address

**OVERVIEW**

All the instruments are used in orthopaedic surgery of the knee or hip. They are for use in an operating room only, under strictly aseptic conditions.

**Stated technical performance**

- Instrument durability: resistance to wear, resistance to mechanical stress when handled by the user under normal conditions, maintaining functionality.
- Optimizing the duration of surgical procedures.
- Biocompatible materials, per NF EN ISO 10993-1.
- The instruments are compatible with the various ranges of implants in both size and shape; the instruments are designed to adapt to the devices they are used to implant.

**Indications**

See the instructions for use of the implant in question.

**Contraindications**

- Allergy to the materials of the implant and/or instrument
- Associated conditions that may compromise the function or implantation of the prosthesis
- Morbid obesity

See the instructions for use of the implant for the specific contraindications to the implant(s) in question.

**Adverse events**

- Pain
- Effusion
- Haemarthrosis
- Infection
- Dislocation
- Trauma or injuries to patients' tissue
- Instrument breakage or wear
- Poorly positioned implant
- Implantation difficult or impossible
- Unequal length of the lower limbs
- Arthrofibrosis
- Baker's cyst
- Allergic reaction
- Metallosis
- Bone fracture
- Implant or instrument falling
- Joint instability
- Decreased range of joint motion
- Prolonged surgical time
- Neurological, vascular, endocrine, metabolic complications
- Loosening
- Postponed procedure
- User injury
- Implant migration
- Surgical revision
- Decreased implant performance
- Decreased bone stock
- Premature breakage or wear of the implant
- Wear to the joint surface

See the instructions for use of the implant for the adverse events directly related to the implant(s) in question.

**CHOICE OF INSTRUMENTS**

To choose the correct instrument, the user should refer to the surgical techniques of the range in question, supplied by DS. The surgeons must know the surgical technique, how to use the instruments and the characteristics of the devices implanted before the surgery is performed.

**DUTIES OF THE HEALTHCARE FACILITY**

DS instruments are supplied in baskets, which are inside aluminium containers, used as transport containers. DS does not guarantee that the containers are effective for ensuring and maintaining instrument sterility. The healthcare facility is responsible for sterilizing the ancillary devices, and keeping them sterile, in compliance with the applicable regulations. The healthcare facility must clean, disinfect and sterilize all instrument kits before returning them to DS. The practitioner should always ensure that the staff, equipment and processes in place to treat DS ancillary devices achieve the expected results. This generally requires regular validation and testing of the processes. Every deviation of the relevant department, compared to the instructions provided, must be thoroughly assessed to determine the effectiveness and potential harmful consequences.

Making any change whatsoever to DS instruments is strictly prohibited. DS, and DS alone, has the expertise to do this. DS accepts no liability in the event of non-compliance with the above.

## **LIST OF MATERIALS THAT THE ANCILLARY DEVICES ARE MADE FROM**

<b>Stainless steel*</b>	<b>Polymers</b>
<b>Alloys</b>	Polyether ether ketone (PEEK)
[Ti6Al 4V]	Polyphenylsulfone (PPSU)
CoCr casting alloy	Polypropylene (PP) H
AL6082 - AL7075	Polyoxymethylene (POM) C
NiT (nitinol)	Silicone (Si)

\* Information on the type of stainless steel can be requested for each instrument that contains it

## **LIMITS AND RESTRICTIONS OF TREATMENT**

Repeated cleaning and sterilization cycles have little impact on the ancillary devices. The end-of-life of an ancillary device is generally determined by wear and damage caused during use, insufficient protection, careless handling or an inappropriate treatment method.

Hydrogen peroxide ( $H_2O_2$ ), nitric acid ( $HNO_3$ ), sulphuric acid ( $H_2SO_4$ ) or oxalic acid ( $C_2H_2O_4$ ) based solutions are not

recommended for treating instruments made from titanium alloys (risk of black discolouration). Similarly, immersing any of our instruments in solutions containing aldehyde or high concentrations of chloride (sodium hypochlorite, etc.) is unauthorized.

## **DATE OF INITIAL MARKET LAUNCH OF THE DEVICES**

- ADES Instruments: 2003 & 2005
- SHIVA Instruments: 2005
- SMS Instruments: 2013
- SYMBOL CUP Instruments: 2014
- SYMBOL STEM Instruments: 2015
- TRICCC Instruments: 1997

## **DATE THAT THE LAST INSTRUCTIONS FOR USE REVISION WAS PUBLISHED**

- October 2019

DEDIENNE Santé thanks you for your trust, and remains available should you require further information.

## **REQUIRED OPERATIONS WHEN THE EQUIPMENT IS RECEIVED**

As the practitioner doesn't know all of the variables involved in handling and transportation, he/she should clean and disinfect the equipment and, if necessary, perform prion inactivation (see the *Prion Inactivation* chapter at the end of this document).

If entering storage, the equipment should be stored in a clean, dry room. The equipment should not come into contact with other devices. The mechanisms and assemblies must be verified: there should not be too much resistance when the instruments are moved, and there should be no marks from damage or excessive wear.

DS instruments are not supplied sterile. It is essential that they are treated (preparation, cleaning, disinfection, drying, sterilization cycle) before use, according to the instructions below.

## **TREATMENT INSTRUCTIONS**

### **Specific instructions after use in the operating room**

The instruments must be treated within two hours of use to minimize the risk that stains dry before cleaning.

Remove excess body fluid and tissue from the instruments, by running them under cold water or using a disposable lint-free wipe.

Place the instruments in a basin of distilled water or, preferably, on a tray covered with moist towels.

To avoid the risk of contamination, used instruments must be transported to the supply centre in closed or covered containers, per best practice for hospital pharmacies.

Before cleaning, don't let saline solution, blood matter, body fluid, bone and tissue fragments or other organic matter dry on the instruments.

### **Preparation before cleaning**

Equipment: detergent/disinfectant (ALKAZYME - Manufacturer: ALKAPHARM - concentration: 0.5%); soft, non-metallic brush; long, narrow, soft-bristled brushes; ultrasound tray; untreated water; distilled water.

Step 1: Disassemble the instruments, if required (see appendix 1)

Step 2: Rinse the instruments under water until macroscopic stains are removed. Pay particular attention to holes and other hard-to-access areas.

Step 3: Immerse the device fully in enzymatic detergent and treat in the US for 15 minutes at 37 kHz.

### **Cleaning, disinfection and drying**

DS has validated two cleaning and disinfection methods, one manual and one automatic. Either may be used, depending on the equipment available in the healthcare facility.

#### ***Manual protocol:***

Step a: Prepare the ALKAZYME detergent/disinfectant solution per the manufacturer's recommendations: immerse the ancillary devices fully, with any hinges open wide. Leave to soak for at least 15 minutes.

Step b: Move mechanisms such as hinges, handles, etc. to remove trapped blood and fragments. Rub the entire surface of

the devices with a soft, non-metallic brush, paying particular attention to grooves and other hard-to-access areas. Rub the inside of tubes and holes with a long, narrow, soft-bristled brush, using a circular motion.

Step c: Rinse the devices with room temperature distilled water for 1 min. 30 sec. Then leave to soak for 1 min. 30 sec., then move the mobile parts, using a syringe for holes and hard-to-access areas. Blind holes should be filled and emptied repeatedly.

Step d: Wipe the instruments with a clean, absorbent, non-woven cloth or clean, lint-free fabric. Dry the inside with medical compressed air.

Step e: Check that the instruments have no visible stains. Repeat the cleaning procedure from step (a) if there are visible stains, and check again.

Check that the instruments have no traces of cleaning product. If there are visible traces, rinse the instrument thoroughly until these traces are removed. Then repeat step (d).

#### **Automatic protocol:**

Step A: Place the instruments in the basket of a suitable washer-disinfector, leaving any hinges open so that water can get into grooves and holes.

Step B: Apply the following cleaning cycle:

	Product type and concentration	Time per step (effective)	Temperature	Water
Steps	Pre-wash	2'	Room temperature	Untreated water
	Wash	MediClean 0.54%	55°C	Untreated water
	Rinse 1	10'	Room temperature	Untreated water
	Rinse 2	2'	Room temperature	Distilled water
	Thermal disinfection	6'	93°C	Distilled water
	Dry	30'	Air: 110°C	
	Cool	5'	Air: 30°C	

Step C: Check that the instruments have no visible stains. Repeat the cleaning procedure from step (A) if there are visible stains, and check again. Check that the instruments have no traces of cleaning product. If there are visible traces, rinse the instrument thoroughly until these traces are removed. Then wipe the instruments with a clean, absorbent, non-woven cloth or clean, lint-free fabric. Dry the inside with medical compressed air.

Once the cleaning/disinfection steps are complete, reassemble the devices.

#### **Tests and maintenance**

Check the mechanisms (sliding pieces, levers, etc.), there should not be too much resistance when they are moved. Examine each device to check that there is no damage and/or excessive wear.

Damaged or misshapen instruments, with signs of shock or oxidation, should be removed and replaced to avoid a risk of injuring the patient. The same applies to sharp edges with nicks.

If polymer surfaces appear "chalky", with superficial deterioration (cracks or delamination, for example) or if polymer devices are misshapen or warped, they should be replaced.

If an instrument should be replaced, report it to your DS representative. All instrument-related incidents must be reported to DS.

#### **Packaging**

To avoid any further contamination, the ancillary devices should be packaged quickly.

- Individual packaging: use an appropriate pouch for steam sterilization. Ensure that it is big enough to hold the ancillary device without pulling the edges. Do not put multiple ancillary devices into the same pouch.
- Batches: ancillary devices should be put in the original basket supplied by DS. Before shipping, these baskets underwent systematic verification, performed by DS (cleaning and visual inspection). Put the ancillary devices back in the correct place, as indicated in the check-list and photos supplied with the containers, to avoid deterioration during transport. Check that the ancillary devices are not touching each other, particularly sharp devices. Pack everything into appropriate packaging for steam sterilization, per best practice for hospital pharmacies.

#### **Sterilization**

The following sterilization method has been validated by DS: steam autoclave at 134°C for at least 18 minutes, then 20 minutes drying time. This process is compliant with ISO 17665-1 "Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices".

When sterilizing several baskets of ancillary devices at once, ensure not to exceed the maximum load of the autoclave.

## **Storage**

Instrument kits that have been treated and packaged should be stored in a cool, dry place.

## **Transport**

To prevent potential damage during transport, using the baskets and containers supplied with the ancillary devices is recommended. Particular attention should be paid to protecting the devices and packaging before shipping.

## **Disposing of devices**

If applicable, the instruments should be disposed of in compliance with the recommendations for eliminating medical waste with an infectious risk.

## **PRION INACTIVATION**

After the cleaning phase, the equipment should undergo prion inactivation, if necessary. The choice of inactivation procedure should be made by the client, per Instruction DGS/R13 no. 2011-449, dated 1st December 2011.

After prion inactivation, the equipment should be systematically rinsed carefully with water. Dry the equipment before packaging.

## **RESHIPPING TO DEDIENNE Santé**

Fill out the FOR\_LOG\_005 treatment traceability sheet accompanying the instruments, then ship the clean, decontaminated equipment to DEDIENNE Santé. **It is mandatory to return the completed traceability sheet.**

## APPENDIX 1: INSTRUMENT DISASSEMBLY PROCEDURE

All instruments provided disassembled in baskets must be disassembled after use, in order to clean and sterilize them. The instruments listed below must be disassembled before cleaning, per the specified instructions, then reassembled and placed in the basket before sterilization.


### 1. SYMBOL CUP Instrum

Name/reference	Disassembly procedure
 Press for mobile liner (3700502204476)	<ul style="list-style-type: none"><li>- Unscrew the head guide pin screw (A) (3700502205817) from the head guide (B) (3700502205800)</li><li>- Unscrew the pusher (C) (3700502205824)</li></ul>
 Bottom impactor tips (3700502204841)	
 Side tip (3700502204827)	
 D30 reduction tip (3700502204803) D40 reduction tip (3700502204810)	Unscrew the plastic part from the metal part of each tip.
 Delta liner gripper tips (3700502205268; 3700502205275; 3700502205282)	
 Suction cup (3700502204834)	




### 2. SYMBOL STEM Instrum




Name/reference	Disassembly procedure
 Gripper tips (3700502205756; 3700502205763; 3700502205770)	Unscrew the metal part from the plastic part of each tip.


### 3. SMS Instrum

Name/reference	Disassembly procedure
 <p>Graft preparation workbench. This consists of:</p> <ul style="list-style-type: none"> <li>- Locking screw (A) (N352062)</li> <li>- Workbench lock nut (B) (N353098)</li> <li>- Graft plate (C) (N353096)</li> <li>- Pressure monitor (D) (N352066)</li> <li>- Mobile ligament tensor (E) (N352060)</li> <li>- Workbench (F) (N352059)</li> </ul>	 <ol style="list-style-type: none"> <li>1. Unscrew the locking screw (A) and the lock nuts (B)</li>    <li>2. Remove the graft plate (C)</li>       <li>3. Remove the mobile ligament tensor (E) and the pressure monitor (D) from the workbench (F)</li> </ol>


### 4. SHIVA Instrum and TRICCC Instrum

Name/reference	Disassembly procedure
 <p>The extractor (A) (0402537), feeder (B) (0402537-2) and slide wheel (C) (0402537-3) should be disassembled before cleaning</p>	 <p>These three instruments are assembled in the basket on arrival. The slide wheel should be unscrewed, then the feeder removed from the extractor.</p>
Tibial cutting block base (96806)	 <p>The tibial cutting block base (A) is assembled with an adjustable tibial cutting lock nut (B) (96927) and screw (C) (96832). These two components need to be unscrewed from the tibial cutting block base during cleaning.</p>
Slider (A) (0304596)	 <p>This instrument is assembled with a wheel (B) (0304597) and a screw (C) (96832), which should be unscrewed.</p>

Wheel (0304570-3)		 Unscrew the wheel from the 4 in 1 femoral cutting block (0304564 to 0304566).
Clamp for patellar cementing (A) (03045107)		 Remove the tightening ring (B) (03045107-10)
Femoral impactor (96834) Tibial impactor (96835) Tibial impactor (0304545)	  	Remove the impactor tips by unscrewing both screws located on the T-shaped section.
Condyle impactor-extractor		Remove the tip (B) of the impactor-extractor by unscrewing the threaded stem (A).
Distal valgus cutting guide (96898)		 Remove the peg (A) from the adjustable holes (B).
Extramedullary sighting tube (A) (96802)		This instrument is assembled with two screws (B) (96832). These screws should be unscrewed before cleaning.

Tibial gauge (A) (0304568) or tibial index (B) (96809-2) Tibial index base (C) (96809-1)	 Disassemble the tibial gauge or index from the tibial index base.
 Safety fork impactor (96836)	Unscrew the yellow tip (A).
 Test condyles (0304526 to 0304528 R & L)	<ul style="list-style-type: none"> <li>- Unscrew the test intercondylar component (A) (0304526-3 to 0304528-3)</li> <li>- Remove the pivot bolt (B) (0304526-2 to 0304528-2)</li> <li>- Remove the locking screw (C)</li> </ul>
 Ankle clamp (96803)	Remove the spring (A) from the bolt (B).

## 5. ADES Instrum

Name/reference	Disassembly procedure
Press (A) (051300), combined with: <ul style="list-style-type: none"> <li>- Press base (B) (051700)</li> <li>- Cone press base (C) (051608; 051610; 051610; 051612)</li> </ul>	 After use, the press base should be removed from the press using a no. 3 Allen key (02011038). Remove the cone press base from its location on the press.
 Cup impactor (050700)	Unscrew the blue tip (A).
 Liner and head impactor (051400)	Unscrew the yellow tip (A).
 Gripper screw (051000)	Unscrew the yellow tip (A) from the gripper screw.

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